

dFlow

Electronic Bottle-Top Dispenser

Features

- Digital operation
- Electronic control provides a fatigue-free operation
- Remote control panel prevents manual disturbance during operation
- Two operation modes:
 - Traditional Dispensing, the liquid is continuously dispensed
 - Stepper Function, for dispensing series



Specifications

Volume Range	0.1-99.9mL, Max piston lift is 10mL, resolution is 100μL
Volume Accuracy	Dispenser: R = 10mL ±30μL; CV = 10ml ±10μL Stepper: R = 1mL ±6μL; CV = 1ML ±9μL
Operating Temperature Range	10°C-30°C
Quality Control	DIN EN ISO 8655
Contro Type	External control